



## TECHNICIAN TALKING POINTS

---

Setting the patient up for a successful AdaptDx Pro<sup>®</sup> experience is simple with this device. In fact, the on-board technician, Theia<sup>™</sup>, does most of the work by coaching the patient through the test. **Below are some examples for how you can address patient questions and set expectations.**

“This test will measure how quickly your vision adapts from light to dark.”

“Here is the controller you will use to respond to the prompts provided by your on-board technician, Theia.”

“You can press any button on the controller to indicate your response.”  
*Hand the controller to the patient.*

“Theia will guide you step-by-step through the testing process. She will give you rest breaks and encouragement to help you stay focused.”

“We need to place the device on your head. I will give you instructions on how to place the headset on your head.”

*Use the following instructions:*

- Ask the patient to place the controller in his/her lap while donning the headset.
- Give the patient the headset, asking them to place one hand on the front and the other on the adjustment dial on the back of the device.
- Make sure the patient can see the red light inside the device.
- Ask the patient to place the eyecups firmly around their eyes, while holding the device in place.
- Slide the back of the headset onto the head and use the dial to adjust the fit until it's secure.
- Engage the top Velcro strap by pulling up and back to ensure the headset is comfortable across the bridge and nose.

“When you are ready to begin the test, press your response button and Theia will take it from there.”

*When the test is complete:* “Would you like my help removing the device?”